

300 mL Reagent Reservoir 96-well, SLAS-ANSI

Specification Sheet (Part/REF # 56669-01)



250 mL SBS Reagent Reservoir 96-well

Product Specifications

Material	Polypropylene - PP
Recycling Code	5
Quantity Per Case	Case of 5 reagent reservoirs
SLAS ANSI	Yes
Barcode	No
Compatible Carrier	SLAS ANSI carrier

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/56669-01>

Spec. Version B