

33 gauge, Small Hub RN Needle, 0.5 in, point style 4 at 12°, 6/PK

Specification Sheet (Part/REF # 7803-15)



Small Hub Removable (RN) needles can be used with Removable Needle (RN) syringes that are 100 µL and smaller. Standard needle length is 2" / 51 mm. Removable (RN) needles are available in convenient six packs. Custom lengths and point styles are available.

Product Specifications

Gauge	33 gauge
Point Style	4: Beveled
Needle Length	0.50 in (12.7 mm)
Autoclavable	Yes
Gas Sterilizable	Yes
Sales Quantity	6 needles per pack
Hub	Small RN
Bevel Length	0.99 mm
Bevel Angle	12°

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/7803-15-2>

Spec. Version B