

300 µL Conductive Sterile Filter Tip Boxes

Specification Sheet (Part/REF # 235261)



Tip Boxes

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS
Tip Type	Standard
Tip Volume	300 µL
Maximum Tip Capacity	300 µL
Filter	Filter
Sterile	Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	Yellow
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Box and Lid Material	Polystyrene
Box Recycling Code	6
Autoclavable	No
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
SLAS ANSI	Yes
Tip Length	59.9 mm
Orifice Size	0.4 mm
Case Size	Case Size of 5,760 tips (6 layers of 10 boxes)
Packaging Option	Individually Wrapped CO-RE Tip Boxes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/disposable-tips/300µl-conductive-sterile-filter-tip-boxes>