

# Semipreparative HPLC Column Repack (7.0, 7.8, 10 mm ID)

Specification Sheet (Part/REF # 79696)



Inherent reversed-phase characteristics make the PRP-1 a rugged, multi-purpose column.

## Product Specifications

<b>Separation Mode</b>	Reversed-Phase Anion Exchange Cation Exchange Ion Exclusion
<b>Stationary Phase</b>	(Choose a phase)
<b>Format</b>	HPLC Column
<b>Hardware Inner Diameter</b>	7.0 mm 7.8 mm 10 mm
<b>Hardware Length</b>	(Choose a length)
<b>Particle Size</b>	(Choose a particle size)
<b>Hardware Material</b>	Stainless Steel

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79642-2>

**Spec. Version B**