

Preparative HPLC Column Repack (100 mm ID)

Specification Sheet (Part/REF # 79699)



Inherent reversed-phase characteristics make the PRP-1 a rugged, multi-purpose column.

Product Specifications

Separation Mode	Reversed-Phase Anion Exchange Cation Exchange Ion Exclusion
Stationary Phase	(Choose a phase)
Format	HPLC Column
Hardware Inner Diameter	100 mm
Hardware Length	(Choose a length)
Particle Size	(Choose a particle size)
Hardware Material	Stainless Steel

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/hplc-columns/79645-2>

Spec. Version B