

Custom HPLC Column Packing Service - Microbore/Analytical (1.0, 2.1, 4.1, 4.6 mm ID)

Specification Sheet (Part/REF # 79651)



Product Specifications

Stationary Phase	Customer Supplied
Format	HPLC Column
Hardware Inner Diameter	Customer Specified 1.0 mm 2.1 mm 4.1 mm 4.6 mm
Hardware Length	Customer Specified 50 mm 100 mm 150 mm 250 mm
Hardware Material	Customer Specified (Stainless Steel or PEEK)

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/hplc-columns/79641-2>

Spec. Version B