

# Intraocular Injector Syringe (6.5 µL), Removable Needle (RN), (Needle Sold Separately)

Specification Sheet (Part/REF # 6609071-01)



Enables one-handed syringe operation for virus or plasmid injections.

## Product Specifications

Needle Sold Separately	Yes
Volume	6.5 µL
Syringe Type	Microliter
Syringe Series	800
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Fluid Path	SST, borosilicate glass
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	2000 psig (137.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.019 inches (0.485 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Std Needle Dead Volume	N/A
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/syringes/7642-01-2>

Spec. Version B