


ZEUS X1 Replacement 1 mL Pipette Head

Specification Sheet (Part/REF # 6606804-11)



ZEUS® X1 air displacement technology employs Hamilton's innovative air displacement pipetting technology and CO-RE II® tips.

Product Specifications

Accuracy	<1% at 1000 µL, <5% at 1 µL
Precision	<0.75% at 1000 µL, <4% at 1 µL
Volume Range	1 - 1,000 µL
Compatible CO-RE®II Tips	10 µL, 50 µL, 300 µL, 1 000 µL

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/oem/6606801-01-2>