

Micro Plate 96-Well for Fluorescence Imaging

Specification Sheet (Part/REF # 10165328)



Hamilton 96-well PCR FramePlate®

Product Specifications

Maximum Well Capacity	392 µL
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Material	Polystyrene
Recycling Code	6
Autoclavable	No
Well Bottom	Flat
SLAS ANSI	Yes
Barcode	Yes
Compatible Carrier	FLUOREYE device or SLAS/ANSI carrier
Quantity	40 plates per case, 10 plates per bag
Shelf Life	4 years after production date
Temperature Range	-20 °C to +60 °C

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/other-robotics/814302-2>

Spec. Version B