

300 µL Nested Conductive Tips

Specification Sheet (Part/REF # 235950)



300 µL Nested Conductive Non-Sterile Non-Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS Microlab Prep
Tip Type	Standard
Tip Volume	300 µL
Maximum Tip Capacity	300 µL
Filter	Non-Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	Yellow
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polystyrene - PS
Rack Recycling Code	6
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
Tip Length	59.9 mm
Orifice Size	0.4 mm
Case Size	Case of 11,520 tips (NTR 5 x 4 stack)
Packaging Option	Nested 96 Tip

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/disposable-tips/300-µl-nested-conductive-tips>

Spec. Version B