

50 µL Clear Tips - Currently unavailable in Europe

Specification Sheet (Part/REF # 235836)



50 µL Clear Non-Sterile Non-Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS Microlab Prep
Tip Type	Standard
Tip Volume	50 µL
Maximum Tip Capacity	50 µL
Filter	Non-Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Clear (Non-conductive)
Barcode Label Color	Light Red
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Tips / Rack	96 Tips/Rack
Tip Length	50.4 mm
Orifice Size	0.4 mm
Case Size	Case of 5,760 tips (Blister 5 x 96 tips per rack)
Packaging Option	Blister 96 Tip

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/disposable-tips/50-ul-clear-tips>

Spec. Version B