

50 µL Clear Filter Tips - Currently unavailable in Europe

Specification Sheet (Part/REF # 235829)



50 µL Clear Non-Sterile Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS Microlab Prep
Tip Type	Standard
Tip Volume	50 µL
Maximum Tip Capacity	50 µL
Filter	Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Clear (Non-conductive)
Barcode Label Color	Light Red
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Tips / Rack	96 Tips/Rack
Tip Length	50.4 mm
Orifice Size	0.4 mm
Case Size	Case of 5,760 tips (Blister 5 x 96 tips per rack)
Packaging Option	Blister 96 Tip

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/disposable-tips/50-ul-clear-filter-tips>

Spec. Version B