

5000 µL Conductive Tips

Specification Sheet (Part/REF # 184020)



5000 µL Conductive Non-Sterile Non-Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS
Tip Type	Standard
Tip Volume	5000 µL
Maximum Tip Capacity	5000 µL
Filter	Non-Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	Green
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	5000 µL channels
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	24 Tips/Rack
Tip Length	118 mm
Orifice Size	1.05 mm
Case Size	Case of 720 tips (Blister 5 x 24 tips per rack)
Packaging Option	Blister 5 mL Tip

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/disposable-tips/5000-µl-conductive-tips-1>

Spec. Version B