

1000 µL Wide Bore 3.2 mm Orifice Conductive Sterile Filter Tips

Specification Sheet (Part/REF # 235842)



1000 µL Wide Bore 3.2 mm Orifice Conductive Sterile Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS
Tip Type	Wide Bore
Tip Volume	1000 µL
Maximum Tip Capacity	1000 µL
Filter	Filter
Sterile	Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	White
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
Tip Length	80.1 mm
Orifice Size	3.2 mm
Case Size	Case of 3,840 tips (Blister 5 x 96 tips per rack)
Packaging Option	Blister 96 Tip

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/disposable-tips/1000-µl-wide-bore-3-2-mm-orifice-conductive-sterile-filter-tips>

Spec. Version B