

1000 µL Wide Bore 3.2 mm Orifice Conductive Sterile Tips - Only available in Europe

Specification Sheet (Part/REF # 235841)



1000 µL Wide Bore 3.2 mm Orifice Conductive Sterile Non-Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS
Tip Type	Wide Bore
Tip Volume	1000 µL
Maximum Tip Capacity	1000 µL
Filter	Non-Filter
Sterile	Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	White
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
Tip Length	80.1 mm
Orifice Size	3.2 mm
Case Size	Case of 3,840 tips (Blister 5 x 96 tips per rack)
Packaging Option	Blister 96 Tip

Specifications are subject to change without notice

Spec. Version B

<https://www.hamiltoncompany.com/disposable-tips/1000-µl-wide-bore-3-2mm-orifice-conductive-sterile-tips>