

1000 µL Wide Bore 3.2 mm Orifice Conductive Filter Tips

Specification Sheet (Part/REF # 235541)



1000 µL Wide Bore 3.2 mm Orifice Conductive Non-Sterile Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS
Tip Type	Wide Bore
Tip Volume	1000 µL
Maximum Tip Capacity	1000 µL
Filter	Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	White
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Tips / Rack	96 Tips/Rack
Tip Length	80.1 mm
Orifice Size	3.2 mm
Case Size	Case of 3,840 tips (Blister 5 x 96 tips per rack)
Packaging Option	Blister 96 Tip

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/disposable-tips/1000-µl-wide-bore-3-2-mm-orifice-conductive-filter-tips>

Spec. Version B