

300 µL Wide Bore 1.55 mm Orifice Conductive Filter Tips

Specification Sheet (Part/REF # 235449)



300 µL Wide Bore 1.55 mm Orifice Conductive Non-Sterile Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS
Tip Type	Wide Bore
Tip Volume	300 µL
Maximum Tip Capacity	300 µL
Filter	Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	Yellow
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
Tip Length	51.5 mm
Orifice Size	1.55 mm
Case Size	Case of 5,760 tips (Blister 5 x 96 tips per rack)
Packaging Option	Blister 96 Tip

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/disposable-tips/300-µl-wide-bore-1-55-mm-orifice-conductive-filter-tips>

Spec. Version B